

### Configuration of Influx

Lasers	Detector	Band Pass Filter	Live/Dead	Fixable	Color
Blue 488nm	A	488/10			SSC
	B	530/40	Sytox Green	Green	FITC/A488/GFP
	C	692/40	7AAD/PI		PerCP-Cy5.5
Green 561nm	A	585/29	Sytox Orange		PE
	B	610/20	7AAD/PI	Red	PE-Texas Red
	C	670/30	7AAD/PI		PE-Cy5
	D	750LP			PE-Cy7
Violet 405nm	A	460/50	DAPI/Sytox Blue	Violet	PacBlue/V450/BV421
	B	520/35	Sytox Blue	Aqua	AmCyan/PacOrange/V500/Q525/BV510
	C	610/20			Q605/BV605
	D	660/20			Q655/BV650
	E	710/50			Q705/BV711
	F	750LP			Q800/BV785
Red 638nm	A	670/30	Sytox Red	Far Red	APC/A633/A647
	B	720/40			A700
	C	750LP		Near IR	APC-Cy7/APC-H7